

Test Ordered	Test	Instructions	Cost US\$
	THYROID ANTIBODY PROFILE (D8T)	If the patient is already on thyroid medication the blood should be drawn 4-6 hours after giving thyroid medication. Collect 3 to 5 mL of blood. Add to the redtop tube or serum separator tube (SST). Allow to clot and centrifuge to separate serum. <b>DO NOT SEND UNSEPARATED BLOOD.</b> Transfer serum to empty redtop or plastic tube. 1.0 mL serum is <b>REQUIRED</b> for the Thyroid Antibody Profile. Test includes T3, T4, free T3, free T4, T3 autoantibodies and T4 autoantibodies.	\$65.00
	TgAA, Add on to D8T	Note: Good for genetic screening. Send extra 0.5 mL separated serum. Note: <b>Allow 3-5 days testing time.</b>	\$27.00
	PROFILE 7200 includes CBC, Differential & Superchem including liver enzymes & D8T	Send 1 mL whole blood in lavender-top tube (EDTA), <b>AND</b> 1.5 mL separated serum as described above for D8T Thyroid Antibody Profile. Include ice pack when shipping, but do not freeze sample.	\$87.00
	D2 includes CBC, Differential & Superchem (SA020)	Send 1 mL whole blood in lavender-top tube (EDTA), <b>AND</b> 1.5 mL separated serum as described above for D8T Thyroid Antibody Profile. Include ice pack when shipping, but do not freeze sample.	\$52.00
	DISTEMPER & PARVO TITERS (SA150)	Send 1.0 mL separated serum.	\$35.00
	RABIES VACCINE TITER (S16685)	Send 1.5 mL serum in red-top or plastic tube. <b>Note: 14 - 21 days testing time.</b>	\$75.00
	VON WILLEBRAND TEST vWD (S17123)	Take 1 mL. separated plasma from a citrate tube (light blue) and transfer plasma to empty redtop or plastic tube. Include ice pack when shipping but do not freeze sample. Clearly Label tube as "PLASMA". <b>Note: 7 days testing time.</b>	\$55.00
	PHENOBARBITAL PROFILE (T750)	Blood sample should be drawn 4-6 hours after giving medication. 1.0 mL serum is required. Please use redtop tube, NOT SST tube.	\$45.00
	BROMIDE PROFILE (T730)	Blood sample should be drawn 4-6 hours after giving medication. 1.0 mL serum is required. Please use red-top tube only, <b>NOT</b> SST tube	\$53.00
	SERUM BILE ACIDS PRE & POST (T220)	1 <sup>st</sup> tube is fasting, 2 <sup>nd</sup> tube is 2-4 hours post-fatty meal. 1.0 mL separated serum in <b>each</b> tube. Label tubes clearly 1 <sup>st</sup> & 2 <sup>nd</sup> .	\$44.00
	URINE BILE ACIDS (85645)	Freshly collected urine, 3 mL.	\$43.00
	PANLEUKOPENIA VACCINE TITER (S16053)	Send 1.0 mL separated serum.	\$35.00

**LABEL TUBES CLEARLY WITH PET'S NAME AND OWNER'S LAST NAME.**

- Ship samples in sturdy cardboard box with plenty of bubble wrap.
- Separated-Serum transferred into a plastic or red top tube is recommended. Send by two day carrier.
- Do not freeze sample. Blue and lavender top tubes require an ice pack with a overnight carrier.
- Send by US Priority Mail, US Express Mail, Fed-Ex or other courier.
- Samples received Monday - Saturday
- HEMOPET will submit samples to ANTECH for analysis and Dr. Dodds will follow with issuing reports & interpretations and recommendations.

Send check or money order along with samples and submission form to:

**HEMOPET, 11330 Markon Drive, Garden Grove, CA 92841**

Credit Cards are acceptable, Please contact Helene Proudfoot (709.745.5373 or [harphill@nl.rogers.com](mailto:harphill@nl.rogers.com)) for details.

Any questions, please call the office - 714-891-2022 Fax: 714-891-2123 Monday thru Friday, 8am - 5pm.